

Compare our products



Voltaren Emulgel P



Voltarol 12Hour Emulgel P 2.32% Gel



Voltarol Pain-Eze Emulgel



Voltarol Emulgel P 180 g



Voltarol 140mg Medical Plaster



Voltarol Heat Patch

FORMULATION

Gel

Gel

Gel

Gel

Plaster

Patch

USE FOR

- For targeted pain relief in osteoarthritis

- For relief of joint pain and swelling (inflammation)

- For targeted relief of muscle pain and swelling (inflammation) including back, neck, shoulder, and for sports injuries

- For targeted pain relief in osteoarthritis

- For relief of muscle pain and swelling (inflammation) including neck and shoulder

- Back and muscle pain relief

BENEFITS

- Pleasant smell, non-sticky: Combines moisturising effect of cream and cooling effect of a gel

- Provides all-day concentrated pain relief; apply twice daily morning & evening

- Up to 3 times more effective pain relief*. Non-sticky: Combines moisturising effect of cream and cooling effect of a gel

- Pleasant smell, non-sticky: Combines moisturising effect of cream and cooling effect of a gel

- Contains an anti-inflammatory, for powerful & targeted deep down relief of pain. Discreet, flexible, thin & odourless. Can be worn under clothing

- Delivers up to 8 hours of continuous pain-relieving heat

DOSAGE

- Amount according to treatment area, apply 3-4 times a day

- Amount according to treatment area. Apply only 2 times a day (morning & evening)

- Amount according to treatment area, apply 3-4 times a day

- Amount according to treatment area, apply 3-4 times a day

- 12 hourly application: Apply once in the morning & evening

- One heat patch per day on the same area of skin

DURATION OF USE

- For up to 14 days

- For up to 14 days

- For up to 7 days

- For up to 14 days

- For up to 7 days

- Single use

ACTIVE INGREDIENT

- Contains diclofenac diethylammonium

- Contains diclofenac diethylammonium

- Contains diclofenac diethylammonium

- Contains diclofenac diethylammonium

- Contains diclofenac sodium.

- Non-medicated

PACKAGE SIZE OPTIONS

- Tubes of 30g, 50g and 100g, the 100g tube also features an easy to open cap designed for osteoarthritis sufferers

- Tubes of 30g, 50g and 100g

- Tubes of 30g, 50g and 100g

- Tubes of 180g

- Pack of 2 and 5 available

- 2 Sachets per Box

*vs a non-medicated gel